EEG Filtering: Techniques and Applications

Callyn Villanueva

I. Introduction

Lectroencephalography (EEG) is a technique used to measure and record the electrical activity of the brain through electrodes placed on the scalp. It is widely used in clinical diagnostics, neuroscience research, and brain-computer interface (BCI) applications. EEG signals are often unfortunately noisy (dare I say, annoyingly noisy), making it essential to filter out unwanted interference to obtain reliable data. Filtering is a crucial step in the preprocessing of EEG signals, enabling accurate interpretation of brain activity by removing noise and artifacts. I've just recently took the initiative of learning the pre-processing factor and oh my god – I have a lot to learn.

II. TYPES OF EEG SIGNALS AND ARTIFACTS

EEG signals are categorized into different frequency bands, each associated with specific brain activities. The main bands include delta, theta, alpha, beta, and gamma waves. Delta waves are typically observed during deep sleep, theta waves in light sleep or relaxation, alpha waves during calm, alert states, beta waves during active thinking, and gamma waves during high cognitive functioning. Understanding these signals is essential for accurate EEG interpretation However, EEG signals are often contaminated by artifacts, which are unwanted signals caused by non-brain activity. These artifacts may arise from eye movements (ocular artifacts), muscle activity (e.g., from facial muscles or eye blinks), electrical interference, and movement-related artifacts. Such noise can obscure the brain's electrical activity, making it necessary to apply effective filtering techniques to ensure the reliability of the data (Jia et al., 2019).

III. HOW EEG FILTERING WORKS: LOW-PASS, HIGH-PASS, BAND-PASS, AND NOTCH FILTERS

EEG filtering techniques play a crucial role in isolating the desired signal from unwanted noise or artifacts. The four primary types of filters used in EEG preprocessing are low-pass, high-pass, band-pass, and notch filters. Each of these filters serves a specific purpose, depending on the nature of the noise and the frequencies of interest. Below is an explanation of how each of these filters works in the context of EEG signal processing.

A. Low Pass Filters

Low-pass filters allow frequencies below a certain cutoff frequency to pass through while attenuating (reducing the amplitude of – I just learned a new word) higher frequencies. This is particularly useful when trying to remove high-frequency noise or artifacts, such as muscle activity (e.g., from facial



Fig. 1. A Typical Recording-to-Publication Data Pipeline, Showing Where Filters Are Applied (https://www.cell.com/neuron/fulltext/S0896-6273




Fig. 2. Both the MNE-Python 0.13 and MNE-C filters have excellent frequency attenuation

muscles) or electrical interference from nearby equipment. These filters are commonly used to eliminate "high-frequency" noise that is not associated with the brain's electrical activity, as EEG signals typically contain frequency components in the range of 0.5 to 100 Hz. For example, in EEG, a low-pass filter with a cutoff frequency of 50 Hz would allow frequencies below 50 Hz to pass while attenuating any signal components above 50 Hz, effectively removing noise from sources like power line interference

B. High Pass Filters

High-pass filters, in contrast to low-pass filters, allow frequencies above a certain cutoff frequency to pass through while attenuating lower frequencies. These filters are typically used to remove low-frequency noise such as baseline drift or slow-moving artifacts. Baseline drift, which can be caused by sweat, electrode impedance changes, or slow movements, is a common problem in EEG recordings. For example, a high-pass filter with a cutoff frequency of 0.5 Hz might be used to remove slow fluctuations or drifts that are not representative of brain activity. By removing these low-frequency signals, high-pass filters can help focus on the brain's more dynamic electrical patterns.

Figure 3.




Fig. 3. Figure 4 The sub-band filtering for raw EEG data. https://www.researchgate.net/figure/The-sub-band-filtering-for-raw-EEG-data_fig2 337644761)

C. Band Pass Filters

Band-pass filters are a combination of low-pass and high-pass filters that allow only a specific range of frequencies (band) to pass through while attenuating frequencies outside of this range. These filters are particularly useful in EEG signal processing because they can isolate specific brain wave bands, such as alpha (8-12 Hz), beta (13-30 Hz), or theta (4-8 Hz), which are associated with particular cognitive or mental states.

For instance, if a researcher is interested in alpha waves, which are typically observed when a person is in a relaxed yet alert state, a band-pass filter with a range of 8-12 Hz could be applied. This would isolate the alpha band while filtering out both lower and higher frequency components, such as delta and gamma waves. Band-pass filters are also useful for removing both high-frequency noise and low-frequency artifacts simultaneously, providing a more precise signal that represents the frequency of interest.

Notch filters are a special type of band-stop filter that are designed to eliminate a very narrow frequency band. These filters are particularly effective at removing specific unwanted frequencies, such as power line noise, which commonly occurs at 50 or 60 Hz depending on the region. For example, in a hospital or laboratory setting, electrical equipment often introduces power line interference at 60 Hz (or 50 Hz in some regions), which can significantly contaminate EEG data. Applying a notch filter with a 60 Hz (or 50 Hz) cutoff can remove this narrowband noise while leaving the rest of the EEG signal intact. Notch filters are essential in environments with high electrical interference and are widely used in clinical EEG applications.

IV. CHALLENGES

EEG filtering is not without challenges. One of the main difficulties is balancing the removal of noise without distorting the true brain signals. Over-filtering can lead to the loss of important information, while under-filtering may allow noise to interfere with signal interpretation. Advanced filtering techniques, such as adaptive filtering or wavelet transforms, are increasingly being employed to address these challenges. Adaptive filtering adjusts filter characteristics based on the signal's properties, making it particularly useful in dynamic environments. Wavelet transforms, on the other hand, break the signal into components at different scales and can more effectively separate signal from noise without causing significant distortion



Fig. 4. Experimental EEG Session https://www.nature.com/articles/s41597-022-01898-y




Fig. 5. EEG Signal - Normal Behavior, Pre-State & Seizure https://asp-eurasipjournals.springeropen.com/articles/10.1186/1687-6180-2014-183

V. REAL WORLD APPLICATIONS

EEG filtering plays a critical role in various clinical and research applications. In clinical settings, clean EEG signals are crucial for the accurate diagnosis and monitoring of neurological conditions, such as epilepsy, sleep disorders, and brain injuries. For instance, filtering out muscle artifacts can enhance the detection of epileptic seizures, which often manifest as sudden changes in brain wave patterns.

In the field of BCIs, filtered EEG signals are essential for decoding brain activity and enabling communication between the brain and external devices. Accurate filtering allows for the reliable control of prosthetics, robotic arms, or computer cursors by interpreting brain signals in real-time. Moreover, filtered EEG is also used to monitor cognitive states, such as focus or relaxation, and to detect abnormalities in mental health conditions like depression or ADHD.

VI. CONCLUSION

Each type of filter plays a critical role in EEG signal processing by targeting and isolating specific frequencies that are either of interest (e.g., brainwave rhythms) or considered noise (e.g., muscle artifacts, electrical interference). By using low-pass, high-pass, band-pass, and notch filters appropriately, researchers and clinicians can enhance the clarity of EEG signals, making them more reliable for diagnostic purposes

and for research into brain activity patterns. These filtering techniques, when applied correctly, help to ensure that the data obtained is both accurate and representative of the underlying brain signals.

REFERENCES

- [1] Alotaiby, T.N., Alshebeili, S.A., Alshawi, T. et al. EEG seizure detection and prediction algorithms: a survey. EURASIP J. Adv. Signal Process. 2014, 183 (2014). https://doi.org/10.1186/1687-6180-2014-183
- [2] Bashivan, P., et al. (2024). A novel approach to EEG signal filtering using adaptive techniques. ScienceDirect. https://www.sciencedirect.com/ science/article/pii/S2666720724001152
- [3] De Cheveigné, A., & et al. (2019). Filters: When, why, and how (not) to use them. Neuron, 102(2), 280-293.
- [4] Jia, H., et al. (2019). Review of filtering methods in EEG signal processing. PMC. https://pmc.ncbi.nlm.nih.gov/articles/PMC6427454/